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Lumi-Glow™ Cell Viability Assay System

(AK1001A, AK1001B, AK1001C, AK1001D)

Description

Lumi-Glow™ cell viability assay system is the homogeneous ready to use reagent for quantitatively detecting ATP levels in the cells. It featured high signal to noise ratio, good reproducibility and stability. The signal half-time is over 2 hours. The homogeneous formulation eliminates the procedures of addition of the reagent after cell lysis, as well as decreases experiment errors generated by multiple additions. Moreover, the product is suitable for high throughput sample measurement due to the stable “glow” type luminescence.

Reagent components

Luciferase, luciferin, buffer were mixed and filled in the 10 ml or 100 ml reagent bottles, three packages are showed below:

| 目录号 | 规格 | 可检测的 96-孔板孔数 | 可检测的 384-孔板孔数 |
|---------|---------|--------------|---------------|
| AK1001A | 10 ml | 100 wells | 400 wells |
| AK1001B | 50ml | 1,000 wells | 4,000 wells |
| AK1001C | 200 ml | 1,000 wells | 4,000 wells |
| AK1001D | 1000 ml | 10,000 wells | 40,000 wells |

Procedures

- Add compounds or controls to the cells, then incubate for appropriate time (the incubation time need be optimized)
- Equilibrate the reagent at RT for 30-60 min till the reagent temperature reaches room temperature
- Equilibrate the cells at RT for 10 min
- Gently mix the reagent through inverting the bottle for several times
- Add 100 µl of reagent to 100 µl of cells in 96-well plate or add 25 µl of reagent to 25 µl of cells in 384-well plate
- Transfer the cells to the corresponding wells in the white or black plate if the above cell plates are transparent. Ignore this step if the cell plates are white/clear bottom or black/clear bottom ones.
- Gently mix the cells
- Read the plate in the appropriate luminescence plate reader after 10-15 min of reagent addition

Storage conditions

Store Lumi-Glow™ cell viability assay kit at -20°C. The product can be stored at RT for 1 week, or 4°C for 8 weeks, or -20 °C for one year, or -80 °C for two years, and can tolerate over 10 times of “freeze and thaw”